

DEPFET for structural biology

WAKATSUKI, Soichi (KEK IMSS, SLAC Photon Science, Stanford Univ.
School of Medicine)

SHIMIZU, Nobutaka (KEK IMSS), MATSUGAKI, Naohiro (KEK IMSS)
USHIRODA, Yutaka (KEK IPNS), TANAKA, Shuji (KEK IPNS)
MIYAKE, Hideki (KEK IPNS), HARA, Koji (KEK IPNS)
HIGUCHI, Takeo (Kavli IPMU)

Grant-in-Aid for Scientific Research (S)

“Construction of novel detector system for structural biology using
ultrafast fine pixel detector”

[Principle Investigator: WAKATSUKI, Soichi]

Needs of new detectors for protein xtallography/solution scattering experiments

Protein Crystallography

Now Efficient data collection

Issue How to measure intensities of very small spots from micro beam/micro xtals

Solution scattering


Now Protein folding & solution dynamics of protein-protein interactions

Issue fastest dynamics: ~msec

Targets: structure and dynamics of large complexes and membrane proteins

- High spatial resolution
- High sensitivity
- Large # of pixels

Required properties of X-ray detectors




- Fast readout
- Large dynamics range



CCD
(Quantum315r)

Detector	Effective pixel size	Readout speed	Dynamic range	Readout scheme
CCD	90 μm	<u>300 msec</u>	<u>10^4</u>	Integration
Hybrid PAD	<u>170 μm</u>	<u>3 msec</u>	10^6	Streamlined





Hybrid PAD
(Pilatus 6M)

Needs for fast & high spatial resolution 2D detector for micro/nano beam


DEPFET (DEpleted P-channel FET)

Detector type	# of pixels	Eff. resolution	Readout speed	Dynamic range	Readout
CCD	16 M	90 μm	300 msec	10^4	Integration
Hybrid PAD	6 M	170 μm	3 msec	10^6	Continuous
DEPFET	8 M	<u>50 $\mu\text{m} \times 55\mu\text{m}$</u>	<u>20 μsec</u> (Continuous)	<u>10^7</u> (Integration mode)	<u>Switchable between integration / continuous</u>

Belle II @ Super KEK B-factory



DEpleted P-channel FET



Upgrade of KEK-B, Nobel Physics Award to Kobayashi & Masukawa (Expected to start operation in 2015)

Major development for structural biology appl.
 ① smaller pixels large area
 ② Fast readout: 20 μs

Sensor chip development by Max Plank Institute, Munchen

①

Prep for developing the DEPFET detector

Preliminary experiments to characterize DEPFET sensors for structural biology applications

① – I High spatial resolution images of diffraction patterns

Check of spatial resolution and peak shape with small beam matched for small xtals


① – II Fast readout data acquisition for solution scattering experiments

Protein folding and photo excitation dynamics followed by time-resolved SAXS with 20 μ sec time resolution




Design optimization of for large-area DEPFET detector system

- Spatial resolution, sensitivity, non-uniformity, dynamic range etc.
- Comparison with commercially available detectors




Small DEPFET sensor
256 × 64 pixels



Pseudo large solid angle data collection using a rotation stage


② Development of large area DEPFET detector

②-I Ultrafast readout system




On-the-fly integration of max 50,000 images (24 bit/pixel)

Software for noise reduction



Crystallographic analysis



Max 1 Gbytes/sec




Fiber optics



Protein dynamics

②-II 8M pixel DEPFET detector based structural/dynamics analysis system

8M pixels with 20 DEPFET sensors



Data acquisition and analysis



③

③ Applications to challenging targets


③-I Structural analysis and dynamics of membrane protein complexes & large complexes

Structure dynamics of complexes involved in photo synthesis and respiratory chain : photo excitation dynamics


Photo excitation




③-II Solution studies of domain association-dissociation dynamics and kinetics of signalosomes



Rahighi et al. *Cell*, 2009




IKK α · IKK β · NEMO complex



Linearly ubiquitylated protein substrate

Ring



20 μ sec time resolution dynamics solution scattering experiments

Structural changes of NEMO in complex with linear ubiquitin chains: NF- κ B signal transduction pathway involved in inflammation, apoptosis, and cancer


Current Studies

Preliminary experiments to characterize DEPFET sensors for structural biology applications

DEPFET early prototype sensor

Pixel size: $20 \mu\text{m} \times 20 \mu\text{m}$

- Sensor thickness: $450 \mu\text{m}$
- Sensor size : $5120 \mu\text{m} \times 1280 \mu\text{m}$
(256×64 pixels)
- Slow readout (~100Hz)







Confirm the response to the X-ray


- Quality of image of diffracted X-ray from Protein
 - X-ray image for the protein crystal
 - X-ray image for the solution scattering
- Position resolution

Experimental apparatus


X-ray beam @ KEK Photon Factory








Scanned with XY-stage to make the image on the detection plane




X-ray Diffraction Image from Protein Crystal

X-ray energy 12.4 keV
(wave length=1.0 Å)
@ BL-5A




Comparison with CCD

Clear advantage!




Comparison with CCD

Clear advantage!



X-ray solution scattering image


- Collagen from Chicken Achilles tendon
 - X-ray energy 8.33 keV (wave length=1.488 Å) @ BL-10C
 - 1-dimentional orientation
 - Lattice spacing: $d=653\text{ Å}$



Detection plane size:
20 x 20 sensors



1d conversion by Circular averaging






n: diffraction order

$$2ds\sin\theta=n\lambda$$

Pixel Positional Response to X-ray photon

- KEK PF BL-5A
- X-ray energy : 12.4 keV, single photon signals
- Pencil beam created by two 10 μm pinhole
(beam size not negligible compared to pixel)
- DEPFET sensor moved by XY-stage
- Count hits of pixels at each position




Smearing in pixel response: $\sim 3\mu\text{m}$
← dominated by the beam size

Sensor parameters and layout


Pixel size: 55um x 55um

Number of rows and columns: 250x768 (same as Belle II)

Thickness: 450um (without thinning)




16 such sensors




Readout

Belle II's configuration

Can be shorter for this application



Further downstream




To be discussed

Work packages

1	Sensor production	HLL		
2	Testing of sensors			
3	ASICs			
4	Assembly			
5	Kapton cable			
6	Kapton attach, wire bond			
7a	DHH/FPGA			
7b	DAQ	KEK/IPMU		
8	Power Supplies			
9	Mechanics and cooling			
10	Project Leader			

Funding

Apr. 2014 – Mar. 2015: 60M¥
Apr. 2015 – Mar. 2016: 60M¥



Money transfer (contract) of >5M¥ requires an intermediate trading company, which takes about 20% of margin.
1€ = 140¥ (which was 100¥ when this project was approved)

Summary

- High spatial resolution and very fast readout is essential for protein structural analyses. DEPFET is one of very promising candidates.
- Basic properties has been successfully tested using a small prototype DEPFET sensor for ILC
 - Clear signal image of diffracted X-ray
 - Good position resolution for single photon
- We need 16 sensors + readout *similar* to Belle II PXD
- We will receive 60M¥+60M¥ in JFY2014 and 15, and can transfer some 285k€+285k€ to DEPFET collaboration
- We need your cooperation to go further